

A Safe Way To Remove Stubborn Stains



Red-Out works on fresh and dried stains without the use of bleach or phosphates.

Red-Out is effective on color residues, dyes of various kinds, wine spots, fruit punches, ink residues, furniture leg spots due to varnish or shellac stains, colored medicines, old coffee spots, and urine stains.

**Kool-Aid®
Furniture Leg Stains**

**Wine
Colored Medicines**

DIRECTIONS FOR USE Shake well before using.
Test fabric in inconspicuous place for color fastness before using.

- 1- Saturate carpet or upholstery fibers with **RED-OUT**.
- 2- Cover spot with a clean towel and apply pressure or non-bleeding weight.
- 3- Inspect spot after 30 minutes. If necessary reapply **RED-OUT** and allow to dry.
- 4- Rinse several times with plain water and absorb into a clean paper towel or towel after each rinse.

See back page for test results



Nilodor, Inc.
10966 Industrial Parkway, N.W. P.O. Box 660
Bolivar, OH 44612-0660 U.S.A.
1-800-443-4321 or (330) 874-1017 FAX (330) 874-3366
www.nilodor.com

RED•OUT Spot & Stain Remover

Red wine test on nylon fiber



A



B



C

A - One direct spray of RED OUT with blotting

B - Two direct sprays of RED OUT with blotting

C - Three direct sprays of RED OUT with blotting

RED•OUT Spot & Stain Remover

Kool-Aid® test on nylon fiber



A



B



C

A - One direct spray of RED OUT with blotting

B - Two direct sprays of RED OUT with blotting

C - Three direct sprays of RED OUT with blotting



Nilodor, Inc.

10966 Industrial Parkway, N.W. P.O. Box 660

Bolivar, OH 44612-0660 U.S.A.

1-800-443-4321 or (330) 874-1017 FAX (330) 874-3366

www.nilodor.com

Quarts - **C327-009** / 6 per case

R.T.U. pH - 8

Odor Control & Facility Care Solutions